



Be Right™



Test strips free chlorine, 0-600 mg/L, 100 tests

Product #:	2890200
ZAR Price:	Contact Hach
Ships within 1 week	

Hach test strips are manufactured with a unique reagent pad at the tip, using proven chemistries based on standard reference methods. Just dip the strip, wait until color development occurs, and compare the color of the reacted pad to the chart printed on the bottle.

Ideal for anyone performing semi-quantitative pre-screenings, spot-checks, or field tests

Simplify testing and get results in less than 1 minute, in the field or lab

Eliminate chemical handling and clean-up

Reduce guesswork by screening samples

***Not for use in sea or brackish water**

Specifications

Footnote:	*mg/L unless otherwise noted; ppb = $\mu\text{g/L}$; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Parameter:	Free Chlorine, High Range
Platform:	Test Strip
Range:	0 - 600 mg/L Cl_2
Smallest Increments Steps:	Steps: 0, 25, 50, 100, 200, 400, 600
Units:	100 tests