



## Ammonia test strips, 0 - 6.0 mg/L, 25 tests

**Product #:** 2755325  
**ZAR Price:** Contact Hach  
**Ships within 2 weeks**

Measuring range: 0 - 0.25 - 0.50 - 1.0 - 3.0 - 6.0 mg/L; test strips are ideal in pre-test screening to detect the presence of materials that might interfere with lab testing.

---

### Specifications

Footnote:	*mg/L unless otherwise noted, ppb = $\mu\text{g/L}$ , ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Number of tests:	25
Parameter:	Ammonia, Nitrogen, low range - As $\text{NH}_3\text{-N}$ , For freshwater
Platform:	Test Strip
Range:	0 - 6 ppm
Smallest Increments Steps:	Steps: 0, 0.25, 0.5, 1.0, 3.0, 6.0 ppm
Test Kit Model:	Test Strips
Units:	25 tests