



Reactive Orthophosphate (as PO₄-) Test Strips, 0-50 mg/L, 50 tests

Product #:

2757150

ZAR Price:

Contact Hach

Ships within 1 week

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Specifications

Parameter:	Phosphorus, Orthophosphate (reactive) - as PO ₄
Platform:	Test Strip
Range:	0 - 50 ppm PO ₄
Smallest Increments Steps:	Steps: 0, 5, 15, 30, 50
Test Kit Model:	Test Strips
Units:	50 tests