



## Bag, sterile, WhirlPak, 207 mL, 100/pk

<b>Product #:</b>	<b>2233199</b>
ZAR Price:	Contact Hach
Ships within 1 week	

### Sampling and Handling

Whirl-Pak Sterile Bag Vol: 207 mL (7oz), Dim (W x L): 9.5 cm x 18 cm, (3.75 x 7") 3.0 mil thick

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection. "Write-on" surface accepts waterproof ink. 24 oz bag perfect for incubating petri dishes in water bath applications.

---

### Specifications

Capacity:	207 mL
Quantity:	100