



Membrane filters millipore S-pak type HC , 47 mm, 0.7 um pore size, pk/150

Product #:

2936200

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Certified for fecal coliform analysis. Larger pore size (0.7 um) and funnel-like pores creates a less taxing environment for stressed organisms: lessens media loss from evaporation and easier filtration of high particulate samples common to effluents. Whi

Specifications

Diameter:	47 mm mm
Material:	Mixed-Cellulose Ester
Pore Size:	0.7 µm
Quantity:	150