



Sample vial, 120 mL w/out thiosulfate, 200/pk

Product #:	2634002
ZAR Price:	Contact Hach
Ships within 3-5 weeks	

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene.

Specifications

Capacity:	Fill line at 100 mL (\pm 2.5 mL)
Diameter:	43 mm
Length:	93 mm
Quantity:	200/pk
Special Features:	Sample Vial, without Thiosulfate