



## Buffered Dilution Water, 99 mL, pk/25

<b>Product #:</b>	<b>1430598</b>
ZAR Price:	Contact Hach
Ships within 3 days	

Magnesium Chloride/Potassium Phosphate Buffer, 99 mL. Buffered Dilution Water (required for testing non-potable water). pk/25.

### Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

### Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

### Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

---

## Specifications

Quantity: /pk