



## Buffered dilution water pillows, pk/25

Product #:	2143166
ZAR Price:	Contact Hach
Available	

**Hach buffered dilution water pillows, when prepared according to the instructions, is formulated to the specifications as described in Standard Methods for the Examination of Water and Wastewater, 9050C.1.a.1 "Buffered water".**

Buffered Dilution Water Pillows, PN 2143166, contains 25 pillows each of magnesium chloride and potassium dihydrogen phosphate. Makes 25 L dilution water (1 liter per 1 magnesium chloride + 1 potassium dihydrogen phosphate pillow).

### **Complete Solutions for Testing Recreational Water**

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

### **Complete Solutions for Testing Wastewater**

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

### **Complete Solutions for Testing Drinking Water**

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

---

## Specifications

Quantity:	25/pk
-----------	-------