



Hach Modified Colitag Sterile 120 mL Sample Bottles, 100/ pk, with Shrink-banded, Polystyrene

Product #:	8888006
ZAR Price:	Contact Hach
Ships within 2 weeks	

Hach Colitag presence-absence tests deliver the fastest EPA approved results in the market

Hach Modified Colitag Sterile 120 mL Sample Bottles with sodium thiosulfate tablet. Can be used with all Hach Modified Colitag 16-hour presence/absence tests. 100/pk. Hach Modified Colitag offers the following benefits:

The Fastest Presence/Absence Test on the Market at 16 Hours

Fastest EPA-approved test for coliforms and *E. coli*. Test as late as 5:00 PM and still have your results by 9:00 AM the next day.

One Product that Supports 16- or 22-hour Test Methods

Allows for a "standard" test that requires less set-up, or a

"rapid" test when necessary. One media for both tests reduces

complexity in certification and inventory.

Results are Valid for 48 Hours

(Please Note: When transported in hot weather, the tab may melt. This does not affect product quality or function).

Saturday's tests can be read on Monday—no need to work all

weekend. Results stay valid for 800% longer than the competition.

Specifications

Material:	Plastic
Quantity:	100 /pk