



Silica Reagent Set for 5500sc

Product #:

6783600

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Complete set of reagents and standards for the 5500sc Silica analyser.

Specifications

Instrument: 5500sc Silica

Parameter: Silica

Storage conditions: Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)