



## pH Electrode for 8362 Ultrapure Water Sensor

**Product #:** 08362=A=2000  
**ZAR Price:** Contact Hach  
**Ships within 1 week**

Part of analogue pH sensor 8362, specially designed for pure water pH measurement.

---

### Specifications

Application:	Ultrapure water
Cable Connection:	S7 connector
Connector:	S7 Connector
Length:	120 mm
Material:	Sensor: Glass
Measuring range:	0 - 14 pH
Reference electrode:	None
Sensor Thread:	Pg 13.5
Special Notes:	pure water pH electrode
Temperature Sensor :	PT100
Warranty:	6 months