



Flow-through cell kit, 230 V

Product #: 4745002

OBSOLETE ITEM
This item is no longer available.

Low-pressure flow-cell kit allows for a constant optical path. Sample flow controlled by a flow valve module.

Obtain colorimetric measurements with Hach spectrophotometers and turbidimeters without individual sample cells. Eliminate errors due to scratched or smudged sample cells. No need for optically matched sample cells. For rapid, convenient processing of large numbers of samples.

Specifications

Compliance Certifications:	TUV/CE
Includes:	Flow-cell stand assembly, flow-cell inlet reservoir, reservoir cover, collection drain assembly, flow-cell assembly, interconnecting tubing, flow-cell light shield
Input Voltage:	230 V AC
Output Voltage:	12V DC at 1 A
Sample Volume:	120 mL purge volume
Weight:	3.2 kg