



## **Glass Rod Secondary Turbidity Standard #0.1 NTU for TU5200, TU5300sc, and TU5400sc Laser Turbidimeters, with RFID**

**Product #:**

**LZY901**

**ZAR Price:**

**Contact Hach**

**Ships within 3-5 weeks**

The glass rod standard for any TU5 series turbidimeter facilitates problem-free routine checks and confirmation of measurements without the need to mix liquid calibration standards.

---

### **Required Accessories**

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)