



Cell Adapter, COD Vial, 16 mm, for DR/2000, DR/2010, DR/3000 Spectrophotometers

Product #:

4479900

 **Hazardous**

OBSOLETE ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Used with Hach spectrophotometers and portable laboratories. Makes it possible to use other cells in instruments which normally take a 25 mL cell. For use with Models DR/2000, DR/2010, DR/3000.