



BODTrak II reagent set w/ bottles

Product #:

2952440

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Includes 6 amber glass bottles, spatula scoop, nutrient buffer solution pillows (pk/50), and potassium hydroxide pellets, 25 g.

Specifications

Application:	BODTrak II
Capacity:	473 mL
Description:	Reagent Set with Bottles
Formula:	Reagent Set with Bottles
Package Type:	Replacement set
Quantity:	(See product description)