



BOD Trak II Manometric BOD apparatus with accessories

Product #:

2952400

ZAR Price:

Contact Hach

Ships within 1 week

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Manometric BOD(n) laboratory measurement station, control unit and stirrer with pressure sensors for 6 bottles.

BOD measurement programmes for 5, 7, 10 days, OECD biodegradation test for 30 days.

Measuring range 0-700 mg/L O₂, variable.

Incl. 6 rubber sleeve stoppers, 6 492-mL brown glass bottles, grease, 6 magnetic stirrers, funnel, 1 pkg LiOH Powder Pillows, 1 pkg nutrient Powder Pillows.

Specifications

BOD Bottle:	473 mL/6 pcs
Capacity:	Six 492 mL bottles
Drift:	< 3 mg/L BOD in 5 days
Input Voltage:	110 - 240 V, 50/60 Hz
Instrument:	BOD Trak II Apparatus
Interface:	RS232
Measurement method:	Manometric
Output Voltage:	24V, UL CSA, and TUV approved
Range:	0 - 700 mg/L
Warranty:	12 months
Weight:	4 kg
What's included?:	BOD Trak II instrument, Power cord, 230 VA harmonized power cord with continental European plug; Seal cups (6x); BODTrack II amber sample bottles (6x); BOD Trak II Magnetic stir bars (6x); Spatula scoop; Nutrient Buffer solution pillows (1 pk); potassium hydroxide pellets (1 container)

What's included?

BOD Trak II instrument, Power cord, 230 VA harmonized power cord with continental European plug; Seal cups (6x); BODTrack II amber sample bottles (6x); BOD Trak II Magnetic stir bars (6x); Spatula scoop; Nutrient Buffer solution pillows (1 pk); potassium hydroxide pellets (1 container)