



## TOC X5 Shaker

**Product #:**

**LQV148.99.00002**

**ZAR Price:**

**Contact Hach**

**Ships within 3-5 weeks**

### Test TOC in your lab

The TOC X5 Shaker is required for the TNT810 and TNT811 TOC TNTplus tests. It purges the Total Inorganic Carbon in the vial, which means you do not have to use or clean any lab glassware.

### Prepare up to 8 TOC TNTplus vials simultaneously

The TOC X5 Shaker allows you to prepare 8 samples in parallel in only 5 minutes, shortening the time spent on your analysis.

**The TOC X5 Shaker is required for the TNT810 and TNT811 TOC TNTplus tests.**

---

## Specifications

Mains Connection 2:	47 Hz - 63 Hz
Mains Connection 3:	30 VA
Manual Languages:	English, French, Spanish, Portugese, Chinese
Parameter:	TOC (Total Organic Carbon)
Power Requirement:	100 V - 240 V Wall Socket Power Adapter
Warranty:	12 months
Weight:	2.9 kg
What's included?:	TOC-X5 Shaker, plug-in power supply unit, non-slip mat

---

## What's included?

TOC-X5 Shaker, plug-in power supply unit, non-slip mat