



## Biogas Starter Kit, H2S04 Full set of reagents, acc. & electrode

Product #:

X91T100

 **Hazardous**

**OBSELETE ITEM**  
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Kit including reagents, titrant, accessories & electrode to run FOS/TAK application with Biogas titrator

Contains Electrode, reagents and accessories for running analysis on biogas fermenters with BIOGAS Potentiometric Titrator.

---

### Specifications

Dimensions (H x W x D):	380 mm x 230 mm x 457 mm
Length:	450 mm
Warranty:	12 months
Weight:	5 kg