



Conversion Kit, Forward Titration Buret Systems to Digital Titrator

Product #: 1930011
ZAR Price: Contact Hach
Ships within 1 week

Forward titration conversion kit with digital titrator. Amperometric titrator and stirrer sold separately.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

Forward or Back Titrations for chlorine available with the amperometric titrator.

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

End-Point easier to determine

The number used for calculation of sample concentration is recorded directly from the Digital Titrator's digital counter. This can be graphed for results for use with the amperometric titrator.

Specifications

What's included?:	Kit includes: digital titrator, PAO standard solution cartridges (2), delivery tubes (pk/5) and manual
-------------------	--

What's included?

Kit includes: digital titrator, PAO standard solution cartridges (2), delivery tubes (pk/5) and manual