



## Reagent Addition Module for SAC850 or SAC950

**Product #:**

**X91T096**



**Hazardous**

**OBSOLETE ITEM**

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

The Reagent Addition Module contains peristaltic pumps used to add reagents as part of a sample preparation process with volume accuracy around 1%. The Reagent Addition Module consists of a peristaltic pump assembly, a reagent addition tip and silicon tubing.

**Minimizes sample preparation time**

**Simplifies the operator's tasks**

**Increases sample throughput**