



Dynamic rinse module for SAC850/SAC950 sample changer

Product #:

X91T107



Hazardous

OBSOLETE ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Specifications

What's included?:

Shower head assembly (p/n X31T166), Pump module assembly (p/n X16T019), Silicone tube (p/n X31T154) & Installation instructions (p/n D26T042)

What's included?

Shower head assembly (p/n X31T166), Pump module assembly (p/n X16T019), Silicone tube (p/n X31T154) & Installation instructions (p/n D26T042)