



PTM450 Photocolorimetric titration LED module w/out power

Product #:

X61T005

ZAR Price:
Available

Contact Hach

Photometric Titration Module, w/o power supply for use with OPT300 Photocolorimetric probe

Photocolorimetric titrations become fast and simple when you add a PTM450 Photocolorimetric Titration Module and OPT300 Photocolorimetric Probe to your TitraLab Titration Workstation.

Make light work of your titrations

Ideal for virtually all colorimetric titrations

Maintenance-free

Easy to handle as any standard electrode

Specifications

Dimensions (H x W x D):	45 mm x 60 mm x 109 mm
Length:	110 mm
Weight:	250 g

Required Accessories

- OPT300 Photocolorimetric probe (Item E71T001)