



Carousel 12 positions, 50 mL & 150 mL beaker, AS1000

Product #: LZE203
ZAR Price: Contact Hach

Carousel tray for 50 mL & 150 mL beakers, 12 positions, for AS1000, without beakers

Specifications

Number of sample positions: 12
Warranty: 12 months

Required Accessories

- Set of beakers, 50 mL, PP (5 pcs) (Item LZE108)
- Set of beakers (5 pcs), 150 mL, PP (Item LZE109)
- Magnetic bar, 20 x 6 mm (5 pcs) (Item LZE136)
- Beaker ring adapter, 50 mL, AS1000 (12 pcs) (Item LZE213)