



Titration Application Pack, Food & Beverage: pH, Total Acidity & SO₂ in wine

Product #: AP0013.AT1222
ZAR Price: Contact Hach
No availability shown

Application Pack for the determination of pH, Total Acidity and Free & Total SO₂ in wine. Suited for Titrilab AT1000 Series, Automatic Titrator Model AT1222.

Specifications

Application:	Food & Beverage Applications
Compatible Instruments:	AT1222
Parameter:	pH, Total Acidity, Free & Total SO ₂ in Wine
Type:	Application Pack
Warranty:	12 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.

What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.