



Titration Application Pack Petrochemistry: TAN (Total Acid Number)

Product #: AP0015.AT1102
ZAR Price: Contact Hach
No availability shown

Application Pack for the determination of the Total Acid Number in Petrochem. Suited for Titrallab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

Specifications

Application:	Petrochemical Applications
Compatible Instruments:	AT1102
	AT1112
	AT1122
	AT1222
Parameter:	TAN (Total Acid Number)
Type:	Application Pack
Warranty:	12 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.

What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.