



## Titration Application Pack Petrochemistry: TBN (Total Base Number)

**Product #:** AP0016.AT1102

**ZAR Price:** Contact Hach

Application Pack for the determination of the Total Base Number in Petrochem. Suited for Titrilab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

---

### Specifications

Application: Petrochemical Applications

Compatible Instruments: AT1102

AT1112

AT1122

AT1222

Parameter: TBN (Total Base Number)

Type: Application Pack

Warranty: 12 months

What's included?: Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.

---

### What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.