



## Titration Application Pack Food & Beverage: Salt

**Product #:** AP0010.AT1112

**ZAR Price:** Contact Hach

No availability shown

Application Pack for the determination of Salt in food products. Suited for Titalab AT1000 Series, Automatic Titrator Models AT1112/AT1122/AT1222.

---

### Specifications

**Application:** Food & Beverage Applications

**Compatible Instruments:** AT1112

AT1122

AT1222

**Parameter:** Salt in food products

**Type:** Application Pack

**Warranty:** 12 months

**What's included?:** Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.

---

### What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.