



Titration Application Pack Biogas: FOS/TAK

Product #: AP0006.AT1102
ZAR Price: Contact Hach

Application Pack for the determination of FOS/TAK (Biogas). Suited for Titrablab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

Specifications

| | |
|-------------------------|---|
| Application: | Environmental Applications |
| Compatible Instruments: | AT1102 |
| | AT1112 |
| | AT1122 |
| | AT1222 |
| Parameter: | FOS/TAC (Biogas) |
| Type: | Application Pack |
| Warranty: | 12 months |
| What's included?: | Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper. |

What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.