



Titration Application Pack Water: pH, Alkalinity & Hardness (ISE)

Product #: AP0005.AT1222
ZAR Price: Contact Hach
No availability shown

Application Pack for the determination of pH, Alkalinity and Calcium & Magnesium Hardness in water. Suited for Titrilab AT1000 Series, Automatic Titrator Model AT1222.

Specifications

Application:	Municipal Applications
Compatible Instruments:	AT1222
Parameter:	pH/Alkalinity & Hardness
Type:	Application Pack
Warranty:	12 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.

What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.