



Titration Application Pack Water: Ca & Mg Hardness (ISE)

Product #: AP0003.AT1112

ZAR Price: Contact Hach

Application Pack for the determination of Calcium and Magnesium Hardness in water. Suited for Titralab AT1000 Series, Automatic Titrator Models AT1112/AT1122/AT1222.

Specifications

Application:	Municipal Applications
Compatible Instruments:	AT1112
	AT1122
	AT1222
Parameter:	Ca & Mg Hardness (ISE)
Type:	Application Pack
Warranty:	12 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.

What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.