



Titration Application Pack Beverage: pH, Total Acidity in Soft Drink/Juice/ Energy Drink Bottlers, Double 25 mL

Product #: AP0026.AT1222

ZAR Price: Contact Hach

Application Pack for the determination of pH, Total Acidity in Soft Drink/Juice/ Energy Drink Bottlers. Suited for Titralab AT1000 Series Automatic Titrator Model AT1222.

Specifications

Application:	Food & Beverage
Compatible Instruments:	AT1222
Parameter:	pH, Total Acidity (Double Syringe)
Syringe:	Double
Type:	Application Pack
Warranty:	12 months
What's included?:	Application pack contains a pH probe, USB key with pre-programmed applications, polypropylene beakers, magnetic stirring bars, bottle stoppers, conical adapter and two (2) high precision glass 25 mL syringes with syringe rings.

The AT1000 Instrument and Reagents are sold separately.

What's included?

Application pack contains a pH probe, USB key with pre-programmed applications, polypropylene beakers, magnetic stirring bars, bottle stoppers, conical adapter and two (2) high precision glass 25 mL syringes with syringe rings. The AT1000 Instrument and Reagents are sold separately.