



Disposable BOD bottle, 300 mL

Product #:	2943101
ZAR Price:	Contact Hach
No availability shown	

Single 300 mL disposable BOD bottle without stopper. Plastic bottle is EPA approved for BOD measurement.

Box of 100 bottles available under article number 2943100. Plastic stoppers (25 pack) available under article number 2943900. Glass stoppers (quantity of 1) available under article number 190901.

Specifications

Body Material:	Plastic
Capacity:	300 mL
Material:	PET
Quantity:	1 each