



## CC70-2 Field Carrying Case & Standards Kit for the IQ340/350

Product #:

CC70-2

 **Hazardous**

**OBSOLETE ITEM**

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

CC70-2 Field Carrying Case & Standards Kit for the IQ340/350. Contains Conductivity Standards Bottles 100mL each of 100 uS, 1000 uS, 10,000 uS, 53 mS and Conductivity Test sample bottle.