



ORP110-GS Standard ORP probe

Product #:

ORP110-GS

OBSOLETE ITEM

This item is no longer available.

Standard ORP electrode, BNC connector, 3' lead, 12 x 110 mm barrel, platinum band

Accurate

For precise electrode standardization, one or two buffers are required to check that the electrode is operating properly.

Easy to use

No calibration required. Standardization used with a pH 7 and pH 4 buffer.

BNC connection

Connects to meters with a BNC connection.

Specifications

Cable Length:	1 m
Dimensions (D x L) :	12 mm x 110 mm
Electrode Type:	Laboratory
Filling Solution:	Non-refillable gel
Length:	110 mm
Material:	Sensor Body: Epoxy
Parameter:	ORP/Redox
Probe Junction :	Ceramic pin
Range:	+/- 1999 mV
Sensor Type :	General Purpose
Technology:	NA
Temperature Range:	-10 - 110 °C
Thermistor :	#
Type:	Benchtop
Waterproof Connector (for Portable Meters):	#