



Backpack for portable meters, large, nylon, wheels, w. cases

Product #:	2990800
ZAR Price:	Contact Hach
Call for ship date	

For comfortably transporting portable meters, accessories, and standard solutions.

Convenient

The answer to hands-free portability of analytical testing equipment, these backpacks combine contemporary styling with innovative design.

Customizable

The adjustable divider system allows you to create customized compartments to organize equipment just the way you like it.

Ergonomic

Lightweight, thick padding material gives superior protection during transportation. A stitched foam back panel makes the backpack well suited for either short or long testing adventures. The sternum buckle system adds support when carrying heavier equipment.

A telescoping handle system and a comfort pad provide additional comfort when used as a backpack. The large backpack can also be used as a wheeled bag or a backpack with in-line skate wheels.

Specifications

Dimensions (H x W x D):	432 mm x 356 mm x 152 mm
Length:	152 mm
Material:	Nylon, plastic textile
Weight:	Capacity: 13.6 kg