



## Backpack for portable meters, small, nylon, with cases

Product #:	2990700
ZAR Price:	Contact Hach
Call for ship date	

For comfortably transporting portable meters, accessories, and standard solutions.

### Convenient

The answer to hands-free portability of analytical testing equipment, these backpacks combine contemporary styling with innovative design.

### Customisable

The adjustable divider system allows you to create customised compartments to organise equipment just the way you like it.

### Ergonomic

Lightweight, thick padding material gives superior protection during transportation. A stitched foam back panel makes the backpack well suited for either short or long testing adventures. The sternum buckle system adds support when carrying heavier equipment.

A telescoping handle system and a comfort pad provide additional comfort when used as a backpack. The large backpack can also be used as a wheeled bag or a backpack with in-line skate wheels.

---

## Specifications

Dimensions (H x W x D):	381 mm x 305 mm x 152 mm
Height:	ID:
Length:	152 mm
Material:	Nylon, plastic textile
Weight:	Capacity: 9.1 kg