



Industrial Starter Kit with Small Backpack

Product #: 2991200

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

The Hach soft-sided backpack boasts an innovative design that securely carries all your testing equipment with you wherever you go.

The perfect solution for water treatment professionals or anyone looking to carry analytical testing equipment without all the hassle!

The answer to hands-free portability for analytical testing equipment, these backpacks combine contemporary styling with innovative design

The answer to hands-free portability of analytical testing equipment, these backpacks combine contemporary styling with innovative design.

Customizable

The adjustable divider systems allow you to create customized compartments to organize equipment just the way you like it.

Ergonomic

Lightweight, thick padding material gives superior protection during transportation. A stitched foam back panel makes the backpack well suited for either short or long testing adventures. The sternum buckle system adds support when carrying heavier equipment.

Specifications

Application:	Industrial Starter Kit
Length:	152 mm
Material:	Nylon, plastic textile
Options:	Small Backpack
Weight:	Capacity: 9.1 kg
What's included?:	Includes: DR 900 Portable Colorimeter, MP-6 Portable Meter (pH, ORP, Conductivity, TDS & Resistivity) Essentials Package, Digital Titrator Universal Kit, Small Backpack with Cases and Reagent Sets for Total Chlorine, Free Chlorine, Nitrite, Iron, Phosphate, Hardness, Alkalinity, Sulfite.

What's included?

Includes: DR 900 Portable Colorimeter, MP-6 Portable Meter (pH, ORP, Conductivity, TDS & Resistivity) Essentials Package, Digital Titrator Universal Kit, Small Backpack with Cases and Reagent Sets for Total Chlorine, Free Chlorine, Nitrite, Iron, Phosphate, Hardness, Alkalinity, Sulfite.