



## Drinking Water Starter Kit with Large Backpack

Product #: 2991000

**OBSOLETE ITEM**  
This item is no longer available.

**The Hach backpack's adjustable divider systems allow ultimate customization for the unique needs of each drinking water tester.**

The perfect solution for field operators or anyone looking to carry analytical testing equipment without all the hassle!

**The adjustable divider system allows you to create customized compartments to organize equipment just the way you like it**  
The answer to hands-free portability of analytical testing equipment, these backpacks combine contemporary styling with innovative design.

### Convenient

The adjustable divider systems allow you to create customized compartments to organize equipment just the way you like it.

### Ergonomic

Lightweight, thick padding material gives superior protection during transportation. A stitched foam back panel makes the backpack well suited for either short or long testing adventures. The sternum buckle system adds support when carrying heavier equipment.

---

## Specifications

Application:	Drinking Water Starter Kit
Length:	152 mm
Material:	Nylon, plastic textile
Options:	Large Backpack
Weight:	Capacity: 13.6 kg
What's included?:	Includes: Pocket Colorimeter II, Chlorine MR/HR System, 2100Q Portable Turbidimeter kit, MP-6 Portable Meter (pH, ORP, Conductivity, TDS & Resistivity) Essentials Package and a Large Backpack with Cases

---

## What's included?

Includes: Pocket Colorimeter II, Chlorine MR/HR System, 2100Q Portable Turbidimeter kit, MP-6 Portable Meter (pH, ORP, Conductivity, TDS & Resistivity) Essentials Package and a Large Backpack with Cases