



M27Ag-9 Metal electrode, Ref. Insert, Ag/AgCl, for ISE25Ca-9, K-9, NH₄ & NO₃ (Radiometer Analytical)

Product #:

E34M001

OBSOLETE ITEM

This item is no longer available.

Metal electrode for Karl Fischer titrations.

The right probe for your application.

To meet your analysis needs, the Radiometer Analytical ion selective electrode range offers high quality sensors for a variety of ionic species in solution.

Specifications

Electrode Type:

M27Ag-9