



CDC241-9 Conductivity cell, 2-pole, non-platinised, screw cap (Radiometer Analytical)

Product #:	E61M008
ZAR Price:	Contact Hach
Available	

Conductivity cell for continuous measurements, titration, viscous media.

The right conductivity cell for your application.

It is important to choose a conductivity cell with the right construction and geometry for your particular application and working conditions. The Radiometer Analytical line offers a variety of conductivity cells to meet a large number of applications.

Specifications

Conductivity Cell:	CDC241-9
Electrode:	2-pole, non-platinised
Special Features:	Conductivity cell for continuous measurements, titration, viscous media.
Warranty:	12 months