



CDC741T Conductivity cell, 2-pole, for sample changer

Product #:	E61M012
ZAR Price:	Contact Hach
Available	

2-pole conductivity electrode with MAB6 plug and glass body. Cell constant 1cm^{-1} . Platinised, with temperature sensor. Max. temperature 100°C , min. immersion depth 10mm. For use with sample changer.

Specifications

Conductivity Cell:	CDC741T
Electrode:	2-pole, platinised
Special Features:	Conductivity electrode with MAB6 plug and glass body. For sample changer.
Warranty:	12 months