



CDC745-9 General-purpose conductivity cell, 2-pole

Product #:	E61M013
ZAR Price:	Contact Hach
Call for ship date	

2-pole conductivity electrode with screw cap and epoxy body. Cell constant 1cm^{-1} . Platinised, without temperature sensor. Max. temperature 100°C , min. immersion depth 14mm.

Specifications

Conductivity Cell:	CDC745-9
Electrode:	2-pole
Special Features:	Conductivity electrode with screw cap and epoxy body.
Warranty:	12 months