



CDC866T General-purpose conductivity cell, 4-pole, Temp.

Product #:	E61M015
ZAR Price:	Contact Hach
Available	

4-pole conductivity electrode with MAB6 plug and epoxy body. Cell constant 1cm^{-1} . Platinised, with temperature sensor. Max. temperature 80°C , min. immersion depth 35mm.

Specifications

Conductivity Cell:	CDC866T
Electrode:	4-pole, platinised
Special Features:	Conductivity electrode with MAB6 plug and epoxy body.
Warranty:	12 months