



PHC2401 Combined pH electrode, Red Rod, ring diaphragm, BNC

Product #:

E16M400

ZAR Price:

Contact Hach

Ships within 1 week

Long-term accuracy and reproducibility with Red Rod technology

General-purpose pH electrode with Red Rod reference system. Measuring range 0 to 12 pH, -10 to 100°C. Liquid junction: annular ring. With BNC plug.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

Reproducible measuring conditions ensured by a special bayonet head

Wide temperature range from -10 to 100°C

Fast and accurate response

Stable signal

Trouble-free liquid junction

Specifications

Electrode Type:	PHC2401
Warranty:	12 months