



M295Ag Metal electrode, silver rod, d=7.5 mm, PTFE body; banana plug (Radiometer Analytical)

Product #:	E34M003
ZAR Price:	Contact Hach
Available	

Metal electrode for silver and halide titration.

The right probe for your application.

To meet your analysis needs, the Radiometer Analytical electrode range features single, double, and combined metal electrodes with platinum, silver, gold, antimony, mercury and glassy carbon sensing elements. This variety of designs can be used for redox measurements, redox titrations, potentiometric techniques, and imposed current potentiometric titration.

Specifications

Electrode Type:	M295Ag
Warranty:	6 months