



PHG301 General-purpose pH electrode, glass

Product #: E11M001
ZAR Price: Contact Hach
Ships within 3-5 weeks

pH glass electrode with Ag/AgCl reference system. Measuring range 0-12 pH, temperature range 0-80°C. Min. immersion depth 12mm.

Specifications

Electrode Type:	PHG301
Length:	103 mm
Material:	Glass
Range:	0 - 12 pH pH
Sensor Type :	Glass
Temperature Range:	0 - 80 °C
Thermistor :	none
Warranty:	6 months
What's included?:	Includes: Radiometer Analytical PHG301-9 pH electrode with Certificate of Conformity.

What's included?

Includes: Radiometer Analytical PHG301-9 pH electrode with Certificate of Conformity.

Required Accessories

- CL114 cable, 1m, for screw cap electrode FX/S7/coax, BNC plug (Item A94L114)
- REF321 Reference electrode, Ag/AgCl (Item E21M002)