



## PHG311-9 Glass pH Electrode for Alkaline Samples, with screw cap, without cable

**Product #:** E11M004  
**ZAR Price:** Contact Hach  
**Ships within 3-5 weeks**

Glass pH electrode with screw cap for alkaline samples. Needs cable.

---

### Specifications

Compatible Instruments:	MeterLab & sensION+ GLP lab
Connector:	Screw cap
Electrode Type:	PHG311-9
Length:	120 mm
Material:	Glass
Parameter:	pH
Range:	0 - 14 pH
Sensor Type :	Glass (for alkaline samples)
Special Features:	Glass membrane for reduced sodium error
Temperature Range:	0 - 80 °C
Thermistor :	none
Type:	Lab, external reference probe needed
Warranty:	6 months
What's included?:	Includes: Radiometer Analytical PHG311-8 pH electrode with Certificate of Conformity.

---

### What's included?

Includes: Radiometer Analytical PHG311-8 pH electrode with Certificate of Conformity.

---

### Required Accessories

- CL114 cable, 1m, for screw cap electrode FX/S7/coax, BNC plug (Item A94L114)
- REF321 Reference electrode, Ag/AgCl (Item E21M002)