



Sension pH combination Gel electrode

Product #: 5193500

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Standard combination pH electrode with a non-refillable gel reference and built-in temperature sensor. Contains fixed cable with 5-pin circular connector.

Specifications

Cable Length:	0.91 m
Connector (text):	5-pin circular Connector
Electrode Type:	Laboratory
Filling Solution:	Non-refillable gel
Length:	152.5 mm
Material:	Sensor Body: Epoxy
Probe Junction :	Glass frit
Range:	0 - 14 pH
Sensor Type :	Glass
Temperature Range:	Continuous use: 0 - 45 °C
Warranty:	6 months