



Sodium ISE analysis package

Product #: 2348100

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Sodium Analysis Package includes Sodium Platinum Series Combination Electrode, two ammonium chloride reference cartridges, 100 and 1000 mg/L standards, reagents, beaker, and manual.

Everything you need for ISE measurements in one package.

Hach's ISE Analysis Packages provide a complete solution for ion selective measurements, including ISE, standard solutions, ionic strength adjustor, and replacement electrolyte cartridges.

Maximum life through reduced maintenance.

Sension Platinum Series probes' replaceable gel reference electrolyte cartridge is a fast and easy maintenance solution that maximizes the life of your electrode.

Accurate

You don't just need an answer, but the right one, and fast.

Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

Specifications

Cable Length:	1 m
Diameter:	12 mm
Electrode Type:	Laboratory combination
Filling Solution:	2595802
ISA Required:	4451569
Length:	240 mm
Material:	Sensor Body: PMP
Probe Junction :	Open
Range:	0.01 - saturation mg/L Na ⁺
Sensor Type :	Glass
Special Features:	Platinum Series
Temperature Range:	0 - 50 °C
Temperature Range 2:	N/A N/A N/A
Thermistor :	-

Warranty: 6 months

Required Accessories

- Ammonium Chloride Reference Cartridge, pk/2 (Item 2595802)
- Sodium/Potassium Ionic Strength Adjustor (ISA) Powder Pillows, pk/100 (Item 4451569)