



Be Right™



Sodium Platinum series combination electrode

Product #: 5192500
ZAR Price: Contact Hach
Available

Sodium Platinum Series Combination Electrode with BNC connector. Electrode includes two electrolyte cartridges and a manual.

Accurate

You don't just need an answer, but the right one, and fast. Optimized stabilization algorithms in Hach meter platforms eliminate operator guesswork and reduce inaccurate measurements. Many electrodes work specifically to deliver accurate results and t

Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

Rugged and Durable

Rely on Hach's durable measuring systems designed to withstand the rigors of outdoor testing and frequent use. Choose from the industry's toughest probes or non-glass electrodes for rugged, waterproof systems that deliver day after day.

Flexible

Need to measure in both the lab and field? Hach offers efficiency with reliable, high-performance systems that easily move with you from the benchtop directly to the testing site. Measure the parameters you need today with the assurance you can expand them tomorrow.

Specifications

Cable Length:	1 m
Diameter:	12 mm
Electrode Type:	Laboratory combination
Filling Solution:	2595802
ISA Required:	4451569
Length:	240 mm
Material:	Sensor Body: TPX
Probe Junction :	Open
Range:	0.01 - saturation mg/L Na ⁺
Sensor Type :	Glass
Special Feature:	Platinum Series
Temperature Range:	0 - 50 °C
Temperature Range 2:	N/A N/A N/A
Thermistor :	-

Required Accessories

- Ammonium Chloride Reference Cartridge, pk/2 (Item 2595802)
- Sodium/Potassium Ionic Strength Adjustor (ISA) Powder Pillows, pk/100 (Item 4451569)